Tryptophan ELISA Kit

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

*For the determination of Tryptophan in human EDTA plasma,
serum and urine*

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1. INTENDED USE	22
2. INTRODUCTION	22
3. PRINCIPLE OF THE TEST	23
4. MATERIAL SUPPLIED	24
5. MATERIAL REQUIRED BUT NOT SUPPLIED	25
6. PREPARATION AND STORAGE OF REAGENTS	25
7. PRECAUTIONS	26
8. SPECIMEN COLLECTION AND PREPARATION	27
9. ASSAY PROCEDURE	27
<i>SAMPLE PREPARATION PROCEDURE</i>	27
<i>TEST PROCEDURE</i>	28
10. EVALUATION OF RESULTS	30
<i>EXPECTED VALUES</i>	31
11. PERFORMANCE CHARACTERISTICS	32
<i>CROSS REACTIVITY</i>	32
<i>PRECISION AND REPRODUCIBILITY</i>	32
<i>SENSITIVITY</i>	32
<i>RECOVERY</i>	33
<i>LINEARITY</i>	33
<i>CORRELATION WITH HPLC-MS</i>	34
12. LIMITATIONS	34
13. REFERENCES	34
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	35

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13. LITERATUR

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydroxylase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *AM J Med Genet B Neuropsychiatr Genet.* 2009 Mar 5; 150B (2): 202-8.

Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. *International Journal of Tryptophan Research* 2009: 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry*. 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophane. *Er Neurol*. 1986 ; 25 (Suppl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab*. 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr*. 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. *J Affect Disord*. 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Relation of serum cholesterol, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. *J Psychiatry Neurosci*. 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; *Neurogastroenterol Motil* 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress*, Volume 11, Issue 3 2008, 198-209.

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Manual

Tryptophan ELISA Kit



*For the determination of Tryptophan in human
EDTA plasma, serum and urine*

For research use only

Valid from 05.12.2012

K 7730

1. INTENDED USE

The Tryptophan ELISA Kit is intended for the quantitative determination of Tryptophan in EDTA plasma, serum and urine. It is for research use only.

2. INTRODUCTION

Tryptophan is an essential amino acid, which can not be synthesized by humans and, therefore, must be obtained from diet. Tryptophan plays an important role in the synthesis of serotonin, as it is the only source for synthesizing serotonin (see figure 1).

Serotonin itself is released from the appendices of a widespread transmitter system with global modulating effects on different areas of life. The activity of this system depends not least on the actual tryptophan level. In some disease patterns a tryptophan deficiency is found, which is induced by an activation of the enzyme indoleamine 2,3-dioxygenase (IDO). The activation of IDO causes shifting the metabolism of tryptophan to the kynurenine pathway and its degradation (see figure 1).


Figure 1: Synthesis of serotonin and melatonin via tryptophan (left) and shifting the metabolism of tryptophan to the kynurenine-pathway via the enzyme indoleamine 2,3-dioxygenase (right)

As a consequence of that, the tryptophan level in plasma is reduced; it is not available for serotonin synthesis and in the following, for the melatonin production. As a result, fluctuations of mood, aggression, sleep and eating disorder or depression can occur.

In recent years, an increasing number of disease patterns has been detected, in which either the tryptophan level is significantly lowered due to the activation of IDO or for other reasons, or which show an adverse ratio of tryptophan-to-kynurenine levels. In the former case, the tryptophan deficiency is associated with an increase in inflammation markers. Some indications are mentioned in the following:

- In case of **irritable bowel syndrome** (IBS) interferone-gamma activates the enzyme IDO, causing a tryptophan deficiency,
- In **obese patients**, inflammation markers, which show a chronic activation of the immune system, are increased as well. IDO is also activated here, and the plasma tryptophan levels are low,
- **Stress** activates IDO likewise and facilitates the degradation of tryptophan *via* the kynurenine pathway.

In addition, gastrointestinal dysfunctions like **fructose intolerance** are associated with difficulties in absorbing tryptophan from the gut.

In all these disease patterns fluctuation of mood, sleep and eating disorder or depression are known.

The symptoms of these diseases can be eased by administration of tryptophan. However, such a treatment requires the knowledge of the tryptophan level, which can be determined appropriately with the now available test.

3. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for tryptophan derivatization. Afterwards, the treated samples and a polyclonal tryptophan-antiserum are incubated in the wells of a microtiter plate coated with tryptophan-derivative (tracer). During the incubation period, the target tryptophan in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The tryptophan in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-bound antibody is inverse proportional to the tryptophan concentration in the sample.

During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-tryptophan antibodies. After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow color is inverse proportional to the tryptophan concentration in the sample; this means, high tryptophan concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. Tryptophan present in the patient samples is determined directly from this curve.

4. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7730MTP	PLATE	One holder with precoated strips	12 x 8 wells
K7730ST	STD	Standards (0, 3, 10, 30, 100, 600 μ M)	6 x 200 μ l
K7730KO1 K7730KO2	CTRL 1 CTRL 2	Controls	2 x 200 μ l
K7730WP	WASHBUF	Wash buffer concentrate (10-fold)	2 x 100 ml
K7730AK	AB	Anti-tryptophan antibody (Iyo.)	2 x 1 vial
K7730K	2.AB	POD antibody (concentrate 200-fold)	65 μ l
K7730CSP	2.ABDIL	Conjugate stabilizing buffer	12 ml
K7730RP	REABUF	Reaction buffer	2 x 100 ml
K7730DR	DER	Derivatization reagent	2 x 13.3 mg
K7730LM	DMSO	Dimethylsulfoxide (DMSO)	2 ml
K7730AP	ASYBUF	Assay buffer concentrate (10-fold)	2 x 3 ml
K7730TMB	SUB	TMB substrate	15 ml
K7730AC	STOP	Stop solution	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 μ m) with an electrical conductivity of 0.055 μ S/cm at 25°C (\leq 18.2 M Ω cm).

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μ l** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The wash buffer concentrate is stable at **2-8°C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month.**
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20°C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use. They can be re-frozen up to 3 times.

- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (13.3 mg) in 0.8 ml DMSO**. Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use**. Discard any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.
- Dilute the **assay buffer concentrate (ASYBUF)** with reaction buffer (REABUF) **1:10** before use: Add **27 ml REABUF** to the **3 ml ASYBUF concentrate** in the bottle, mix well. The ELISA kit can be separated into two performances by providing 2 x 3 ml assay buffer concentrate. Diluted assay buffer can be stored at **2-8°C for 4 weeks**.
- Dissolve the content of one vial of **anti-tryptophan antibody (AB)** in **5.5 ml of diluted wash buffer**. The ELISA kit can be separated into two performances by providing two AB vials. Diluted anti-tryptophan antibody can be stored at **2-8°C for 4 weeks**.
- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL) (e.g. **60 µl 2.AB + 12 ml 2.ABDIL, prepare only the required amount**). The undiluted POD antibody (2.AB) is stable at 2-8°C until the expiry date stated on the label. Diluted POD antibody can be stored at **2-8°C for 5 days**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at 2-8°C.

7. PRECAUTIONS

- For research use only.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA plasma, serum and urine

- Serum, EDTA-Plasma and urine are suited for this test system. The samples are stable for up to 6 hours at 2-8°C. For longer storage, blood- and urine samples must be frozen at -20°C. We recommend acidifying the urine samples.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- Samples with visible amounts of precipitates should be centrifuged.
- The EDTA plasma, serum and urine samples are diluted for derivatization (see sample preparation procedure).
- For sample preparation, a derivatization reagent (DER) for derivatization of tryptophan is added (details are given in the sample preparation procedure).

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

Sample preparation procedure

Dilute standars (STD), controls (CTRL) and samples (SAMPLE) by factor **1:41** as follows:

25 µl STD/CTRL/SAMPLE + 1000 µl REABUF (reaction buffer)

Please note: Samples from patients with tryptophan supplementation (e.g. in depletion studies) probably require further dilution, we recommend diluting these samples additionally 1:2 with REABUF.

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials) as follows:

1.	Bring all reagents and samples to room temperature (18-26°C) .
2.	Add 100 µl of diluted standards/ controls/ samples (STD/ CTRL/ SAMPLE) in the corresponding vial.
3.	Add 25 µl of freshly prepared derivatization reagent (DER) into each vial (STD, CTRL, SAMPLE), mix well and incubate for 45 min at room temperature (18-26°C) on a shaker (180-240 rpm).
4.	Afterwards add 1000 µl of assay buffer (ASYBUF) into each vial, mix well and incubate for 45 min at room temperature (18-26°C) on a shaker (180-240 rpm).

2 x 50 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5.	Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) in duplicate on a protocol sheet.
6.	Take as many microtiter plate strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
7.	Wash the microtiter plate strips 5 times by dispensing 250 µl of diluted wash buffer (WASHBUF) into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.
8.	For the analysis in duplicate, take 2 x 50 µl of derivatized standards (STD)/ controls (CTRL)/ samples (SAMPLE) out of the vials and add into the respective wells of the microtiter plate.
9.	Add 50 µl of diluted anti-tryptophan antibody (AB) into each well. Cover the plate tightly.

10.	Incubate overnight (15-20 hours) at 2-8°C .
11.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.
12.	Add 100 µl of diluted POD antibody (2. AB) into each well.
13.	Cover plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal shaker (180-240 rpm).
14.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.
15.	Add 100 µl of TMB substrate (SUB) into each well.
16.	Incubate for 8-12 min at room temperature (18-26°C) in the dark*.
17.	Add 100 µl of stop solution (STOP) into each well, mix thoroughly.
18.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

10. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls, and samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results.**

Exception: The results from additionally diluted samples must be multiplied by this dilution factor.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. *4-parameter-algorithm*

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

2. *Point-to-point-calculation*

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. *Spline-algorithm*

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).


Plausibility of the measured pairs of values should be examined before automatically evaluating the results. If this option is not available in the used program, the pairs of values should be controlled manually.

Controls

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve in $\mu\text{mol/l}$. In the following an example of a calibration curve is given, do not use it for the calculation of your results.

Example of calibration curve

Expected values

Based on internal studies with plasma samples of evidently healthy persons ($n=40$) a mean value of $54 \mu\text{mol/l}$ was estimated. The standard deviation was $12 \mu\text{mol/l}$.

EDTA plasma mean value $\pm 2x$ standard deviation: **$54 \pm 24 \mu\text{mol/l}$**

Normal range: **$30 - 78 \mu\text{mol/l}$**

We recommend each laboratory to develop its own normal range. The values mentioned above are indicative only and can deviate from other published data.

11. PERFORMANCE CHARACTERISTICS

Cross reactivity

5-HTP (5-Hydroxy-Tryptophan)	< 0.5 %
L-Phenylalanine	< 0.1 %
L-Tyrosine	< 0.1 %

Precision and reproducibility

Intra-assay (n=12)		
Sample	Tryptophan [$\mu\text{mol/l}$]	coefficient of variation (CV) [%]
1	17.7	11.4
2	56.7	4.1

Inter-assay (n=5)		
Sample	Tryptophan [$\mu\text{mol/l}$]	coefficient of variation (CV) [%]
1	19.4	6.9
2	59.1	4.0

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 20 times.

Sample	Tryptophan mean value [OD]	2 x standard deviation (SD)	Detection limit [$\mu\text{mol/l}$]
Zero-standard	2.3	0.05	7

Recovery

One serum sample was spiked with different tryptophan concentrations and measured in this assay. The analytical recovery rate was determined by the expected and measured tryptophan levels. The mean recovery rate for all concentrations was 100.8 % (n=5).

Spike [μmol/l]	Tryptophan expected [μmol/l]	Tryptophan measured [μmol/l]	Recovery [%]
0		44.8	
25	69.8	74.7	107.1
50	94.8	89.5	94.4


Linearity

The linearity of the ELISA was determined by the dilution of a spiked serum sample. The mean linearity was 99.8 % (n=5).

Dilution	Expected [μmol/l]	Measured [μmol/l]	Recovery [%]
original		44.8	
1+1	22.4	22.5	100.4
1+3	11.2	11.1	99.1

Correlation with HPLC-MS

10 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.98$.

12. LIMITATIONS

Samples with phenylalanine concentrations above 600 µmol/l should be further diluted and re-assayed.

Strong hemolytic and lipemic samples often show wrong concentrations. Do not measure hemolytic and lipemic samples.

13. REFERENCES

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphysmus is assoziated with diagnosis, but not suicidal behavior, in borderline personality disorder. AM J Med Genet B Neuropsychiatr Genet. 2009 Mar 5; 150B (2): 202-8.

Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. International Journal of Tryptophan Research 2009: 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry*. 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophane. *Er Neurol*. 1986 ; 25 (Suppl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab*. 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr*. 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. *J Affect Disord*. 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Relation of serum cholesterol, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. *J Psychiatry Neurosci*. 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; *Neurogastroenterol Motil* 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress*, Volume 11, Issue 3 2008, 198-209.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number